

Capre F Treatment table

Versatility for physiotherapists

The modifiability, adjustability and silent operation of the sturdy Lojer Capre F improve functionality and enhance the comfort of use. Its multiple height settings, easy-to-use armrests and side supports, the stepless adjustable extension and flexion angle of the foot section improve the ergonomics and efficiency of the work being carried out. All Lojer Capre treatment tables come with high-quality, durable SKAI upholstery.

Lojer Capre – Leading the Way in Treatment Table Safety and Design



Lojer Capre F3
with optional accessories



Lojer Capre F5
with optional accessories



Check out the new Capre product video: <http://bit.ly/capre2nd>

LOJER® For easy care

| Lojer Capre F | Capre F3 | Capre F3H | Capre F5 | Capre F5H |
|---|-----------|-----------|-----------|-----------|
| Table sections | 3 | 3 | 5 | 5 |
| Width (cm) | 55,65,75 | 55,65,75 | 55,65,75 | 55,65,75 |
| Upholstery SKAI Fashion (16 options) | • | • | • | • |
| Upholstery SKAI Clinica (3 options) | option | option | option | option |
| Anatomic design upholstery | option | option | option | option |
| Fully adjustable concave arm rests (gas springs) | - | - | • | • |
| Turning arm rest (can be turned to side supports) | option | option | option | option |
| Side supports with fast locking and horizontal adjustment | option | option | option | option |
| Head support adjustment with gas spring | • | • | • | • |
| Head section, oval shape | option | option | option | option |
| Head section, U-shape | option | option | option | option |
| Double hinged head support | - | - | option | option |
| Leg raise with gas spring | +70°/-12° | +70°/-12° | +70°/-12° | +70°/-12° |
| Drainage with gas spring | option | option | option | option |
| Drainage with electric adjustment | option | - | option | - |
| Central locking castors (Ø75mm) | option | option | option | option |
| Foot adjustment (levelling) | • | • | • | • |
| Height adjustment range (cm) | 45-92 | 49-90 | 45-92 | 49-90 |
| Electric height adjustment with foot pedal | • | - | • | - |
| Hydraulic height adjustment with foot pedal | - | • | - | • |
| Additional hand control or foot pedal | option | - | option | - |
| Height adjustment with 360° foot bar | option | - | option | - |
| Paper roll holder (50 or 60 cm) | option | option | option | option |
| Battery operation | option | - | option | - |
| Safe working load (SWL) | 210 kg | 210 kg | 210 kg | 210 kg |
| CE-marking according to 93/42/EEC | • | • | • | • |
| Warranty (years) | 2/10 y | 2/10 y | 2/10 y | 2/10 y |
| Country of Manufacture | Finland | Finland | Finland | Finland |

Explanation of markings: • =standard feature option =optional accessory - =not available



Capre F5 fully adjustable arm rests



Check out the new feature videos of Capre:
<http://bit.ly/caprefeatures>



Capre F3/F5 the stepless adjustable extension and flexion angle of the footsection



Drainage with gas spring or electric motor



Side supports with fast locking and horizontal adjustment



Upholstery options F3, F5, F5 hourglass



See our upholstery options

<http://www.lojer.com/service/material-options>

LOJER® For easy care

Lojer Oy, Putajantie 42, P.O. Box 54, FI-38210 Sastamala, Finland
 Tel. +358 (0)10 830 6700 | info@lojer.com | www.lojer.com